

Test Inks ORGANIC in bottles

from 31 to 46 mN/m // not subject to labelling
with accuracy measurement $\pm 0,5$ mN/m

STANDARD: 32 34 36 38 40 42 44 46 mN/m
SPECIAL (SE): 31 33 35 37 39 41 43 45 mN/m

STANDARD Set A 32 34 36 38 40 42 44 mN/m 7 glass bottles of 10 ml art.no. 40.20000.4 without brush	OPTIONAL 7 Set 31 to 46 mN/m 7 glass bottles of 10 ml of your choice art.no. 40.20001.4 without brush	PET 10 ml 31 to 46 mN/m especially for use in the food sector art.no. 40.204XX.4 without brush
10 ml 31 to 46 mN/m standard and special art.no. 40.201XX.4 without brush	100 ml 31 to 46 mN/m standard and special art.no. 40.202XX.0	250 ml 31 to 46 mN/m standard and special art.no. 40.203XX.0

Test Pens ORGANIC

from 31 to 46 mN/m // not subject to labelling
with accuracy measurement $\pm 1,0$ mN/m

STANDARD: 32 34 36 38 40 42 44 46 mN/m
SPECIAL (SE): 31 33 35 37 39 41 43 45 mN/m

STANDARD Set 32 34 36 38 40 42 44 46 mN/m 8 test pens art.no. 40.25000.0	OPTIONAL Set of 8 31 to 60 mN/m 8 test pens of your choice art.no. 40.25001.0	OPTIONAL Set of 6 31 to 46 mN/m 6 test pens of your choice art.no. 40.25002.0	OPTIONAL Set of 4 31 to 46 mN/m 4 test pens of your choice art.no. 40.25003.0	Test Pen 31 to 46 mN/m standard and special art.no. 40.451XX.0
---	--	--	--	---

Test Pens ORGANIC Jumbo

32 34 36 38 40 42 44 46 mN/m
with accuracy measurement $\pm 1,0$ mN/m // not subject to labelling

STANDARD Set 30 32 34 36 38 40 42 44 mN/m 8 test pens art.no. 40.26000.0	OPTIONAL Set of 8 30 to 44 mN/m 8 test pens of your choice art.no. 40.26001.0	OPTIONAL Set of 6 30 to 44 mN/m 6 test pens of your choice art.no. 40.26002.0	OPTIONAL Set of 4 30 to 44 mN/m 4 test pens of your choice art.no. 40.26003.0	Test Pen Jumbo 30 to 44 mN/m standard and special art.no. 40.261XX.0
---	--	--	--	---

Test Inks BLUE in bottles

from 18 to 105 mN/m // toxic from 24 to 57 mN/m //
with accuracy measurement $\pm 0,5$ mN/m
30 to 72 mN/m according to ISO 8296 (DIN 53364 and ASTM)

STANDARD: 28 32 35 38 41 44 48 56 mN/m
SPECIAL (SE): 18 (colourless) 20 22 24 26 29 30 31 33 34 36 37 39
40 42 43 45 46 47 49 50 51 52 53 54 55 57 58
60 62 64 66 68 70 72 mN/m (colourless:); 76 84 90 105 mN/m

BLUE Set A 28 35 38 41 44 48 56 mN/m 7 glass bottles of 10 ml art.no. 40.30001.0 with brush art.no. 40.30001.4 without brush	BLUE Set B 28 32 35 38 41 44 48 mN/m 7 glass bottles of 10 ml art.no. 40.30000.0 with brush art.no. 40.30000.4 without brush	BLUE Set C 30 32 34 36 38 40 42 mN/m 7 glass bottles of 10 ml art.no. 40.30003.0 with brush art.no. 40.30003.4 without brush	BLUE OPTIONAL Set 18 (colourless) 20 to 72 mN/m 7 glass bottles 10 ml of your choice art.no. 40.30002.0 with brush art.no. 40.30002.4 without brush
BLUE 10 ml 18 (colourless) to 72 mN/m standard and special art.no. 40.301XX.0 with brush art.no. 40.301XX.4 without brush	BLUE 100 ml 18 (colourless) to 72 mN/m standard and special art.no. 40.302XX.0 art.no. 40.302XX.0 (SE)	BLUE 250 ml 18 (colourless) to 72 mN/m standard and special art.no. 40.303XX.0 art.no. 40.303XX.0 (SE)	COLOURLESS 10 ml 76 84 90 105 mN/m special art.no. 40.301XX.0 (SE) with brush

Test Pens BLUE

from 28 - 60 mN/m in steps of one by one and from 62 - 72 mN/m in steps of two by two
toxic from 24 bis 57 mN/m // with accuracy measurement ± 1 mN/m

BLUE OPTIONAL Set of 8 28 to 72 mN/m 8 test pens of your choice art.no. 40.35001.0	BLUE OPTIONAL Set of 6 28 to 72 mN/m 6 test pens of your choice art.no. 40.35002.0	BLUE OPTIONAL Set of 4 28 to 72 mN/m 4 test pens of your choice art.no. 40.35003.0	Test Pens BLUE 28 to 72 mN/m art.no. 40.351XX.0
---	---	---	---

Test Inks PINK in bottles

from 22 to 60 mN/m // non-toxic STANDARD: 30 32 34 36 38 40 42 44 46 48 50 52 54 56 58 60 mN/m
with accuracy measurement $\pm 0,5$ mN/m SPECIAL (SE): 22 24 26 28 31 33 35 37 39 41 43 45 47 49 51 53 55 57 59 mN/m

STANDARD Set A

32 34 36 38 40 42 44 mN/m
7 glass bottles of 10 ml
art.no. 40.60000.0 with brush
art.no. 40.60000.4 without brush

OPTIONAL 7 Set

22 to 60 mN/m
7 glass bottles of 10 ml of your choice
art.no. 40.60001.0 with brush
art.no. 40.60001.4 without brush

PET 10 ml

32 to 44 mN/m
especially for use in the food sector
art.no.40.700XX.0 with brush
art.no. 40.700XX.4 without brush

10 ml

22 to 60 mN/m
standard and special
art.no. 40.601XX.0 with brush
art.no. 40.601XX.4 without brush

100 ml

22 to 60 mN/m
standard and special
art.no. 40.602XX.0
art.no. 40.602XX.0 (SE)

250 ml

22 to 60 mN/m
standard and special
art.no. 40.603XX.0
art.no. 40.603XX.0 (SE)

Test Pens PINK

from 22 to 60 mN/m // non-toxic STANDARD: 30 32 34 36 38 40 42 44 mN/m
with accuracy measurement $\pm 1,0$ mN/m SPECIAL (SE): 22 24 26 28 31 33 35 37 39 41 43 45 - 60 mN/m

STANDARD Set

30 32 34 36 38 40 42 44 mN/m
8 test pens
art.no. 40.45001.0

OPTIONAL Set of 8

22 to 60 mN/m
8 test pens of your choice
art.no. 40.45000.0

OPTIONAL Set of 6

22 to 60 mN/m
6 test pens of your choice
art.no. 40.45002.0

OPTIONAL Set of 4

22 to 60 mN/m
4 test pens of your choice
art.no. 40.45003.0

Test Pen

22 to 60 mN/m
standard and special
art.no. 40.451XX.0

Test Pens PINK Jumbo

30 32 34 36 38 40 42 44 mN/m // non-toxic
with accuracy measurement $\pm 1,0$ mN/m

STANDARD Set

30 32 34 36 38 40 42 44 mN/m
8 test pens
art.no. 40.46000.0

OPTIONAL Set of 8

30 to 44 mN/m
8 test pens of your choice
art.no. 40.46001.0

OPTIONAL Set of 6

30 to 44 mN/m
6 test pens of your choice
art.no. 40.46002.0

OPTIONAL Set of 4

30 to 44 mN/m
4 test pens of your choice
art.no. 40.46003.0

Test Pen Jumbo

30 to 44 mN/m
standard and special
art.no. 40.461XX.0

QUICKTEST 38®

QUICKTEST 38®

38 mN/m approx.
quick check for polyolefins
5 ml
art.no. 40.55100.0

QUICKTEST 38® JUMBO

38 mN/m approx.
quick check for polyolefins
15 ml, 15 mm line width
art.no. 40.55100.4

RAPIDTEST 38®

RAPIDTEST 38®

38 mN/m approx.
quick check for polyolefins
5 ml
art.no. 40.66100.0

RAPIDTEST 38® JUMBO

38 mN/m approx.
quick check for polyolefins
15 ml, 15 mm line width
art.no. 40.66100.4

EQUIPMENT

COTTON TIPPED APPLICATORS

approved, 100 piece/
pack, 6 inch long
art.no. 40.31700.0

TEST LIGHT

for use in case of low
contrast
art.no. 40.31600.0

arctestCLEANER

Optimizes cleanliness
Content 250 ml
art.no. 40.32000.0

TEST INK CASE SMALL

for 7 test ink bottles of 10 ml, with
storage for test pen set of 4 or 1
pack of cotton tipped applicators
without contents
art.no. 40.31800.0

TEST INK CASE BIG

for 24 test ink bottles of 10 ml
without contents
art.no. 40.31900.0